# STAINS FOR RESIN TEETH

INSTRUCTIONS FOR USE

---

GEBRAUCHSINFORMATION

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## INTRODUCTION

The Stains for Resin Teeth were created to add individual touches to the highly esthetic areas of tooth restorations. To assist you as best possible in your endeavors in the laboratory, CANDULOR offers you a comprehensive solution with this special composition of bonding agent, intense color stains and abrasion-resistant coating layer. The level of individual freedom can be applied to every prefabricated tooth using the Stains for Resin Teeth. The light-curing composite material offers excellent physical and esthetic properties for this purpose.

## PRODUCT DESCRIPTION

The Stains for Resin Teeth Set includes everything you need to convert prefabricated resin and composite teeth into individual restorations easily and quickly. To obtain perfect bonding, the characterized surface needs to be treated with the Composite included in the set. Then the esthetic touches can be added using the different intense color stains. It is also possible to make corrections to the shape using the transparent Finisher paste.

Finally, the characterized surface is covered evenly with the Finisher to ensure best possible protection of the characterized surface against chemical and physical influences.

### Indication

- Characterization of resin teeth using intense color stains
- Protection of the characterized surface with Stains for Resin Teeth finisher
- Shape corrections with the Stains for Resin Teeth finisher

### Contraindication

- Shape modification of anterior incisal edges with pronounced protrusional load
- Characterization of ceramic teeth

## EN

### Advantages

- Natural fluorescent behavior
- Easy adaption through proven consistency of the paste
- Ready-mixed
- Easy handling and targeted application via presentation in a syringe with corresponding tip
- Wide range of hues in shade possible by mixing individual colors

PROPERTIES	STAINS FOR RESIN TEETH*
Modulus of elasticity [MPa]	6500 ± 500
Flexural strength [MPa]	90 ± 10
Hardness (Vickers) [MPa]	440 ± 10
Water absorption [µm/mm <sup>3</sup> ]	15 ± 1
Water solubility [µm/mm <sup>3</sup> ]	1 ± 0

\*The values correspond to examples of test results and do not represent a product specification.

### Warnings

This material has been developed for dental use and must be processed as directed in the instructions for use. The manufacturer does not accept liability for damages caused by use other than directed or improper processing. The user is also responsible for confirming before use that the material is suitable for the intended purpose, particularly if the purpose is not listed in the instructions for use. Avoid contact of uncured Stains for Resin Teeth with skin/mucosa and eyes. Stains for Resin Teeth can be slightly irritating in uncured condition and lead to sensitization against methacrylates. Do not inhale grinding dust. Conventional medical gloves do not offer protection against the sensitizing effect of methacrylates.

**Storage instructions**

- Seal syringes immediately after use (light leads to premature polymerization)
- After finishing remove hollow needle and attach screw cap
- Observe instructions for storage and expiry date on the primary packaging
- Store at temperature 2–28 °C / 36–82 °F
- Protect from sunlight and heat!
- Store out of reach of children

**Processing time**

At a light intensity of 3000 Lux the processing time is approx. 2–25 minutes. The time also depends on the color used. Bright colors react faster than dark ones.

**Delivery form**

- 1 × Stains for Resin Teeth composite, 2 ml
- 1 × Stains for Resin Teeth clear, 1 ml
- 10 × Stains for Resin Teeth intensive colors, 1 ml
- 1 × Stains for Resin Teeth finisher, 2.5 g
- 1 × Instruction for use
- 1 × Tips (contents 20 pieces)

*Delivery form can vary by country.*

**Composition**

Stains for Resin Teeth intense colors consist of:

- > Dimethacrylate (47–48 wt. %)
- > Copolymer and silicon dioxide (49–50 wt. %)
- > Also contains stabilizers, catalysts and pigments (2–3 wt. %)

Stains for Resin Teeth composite consists of:

- > Dimethacrylate (81–82 wt. %)
- > Silicon dioxide (16–17 wt. %)
- > Also contains stabilizers, catalysts and pigments (< 2.5 wt. %)

Stains for Resin Teeth finisher consists of:

- > Dimethacrylate (17–19 wt. %)
- > Copolymer and silicon dioxide, of which inorganic fillers (82–83 wt. %)
- > Catalysts and pigments (64–65 wt. %)
- > Catalysts and pigments (< 1 wt. %)

**POLYMERIZATION PARAMETERS**

Device	Manufacturer	Stains	Composive/ Finisher	Final poly- merization
Quick Lumamat 100	Ivoclar Vivadent AG	20 s Quick	20 s Quick	P2/11 min
Spectramat	Ivoclar Vivadent AG	2 min	5 min	5 min
Laborlight LV-III	GC	2 min	2 min	5 min
Solidilite V	Shofu	1 min	1 min	5 min
Visio beta	3M	4 × 20 s Visio Alfa	4 × 20 s Visio Alfa	2 × 7 min without vakuum
HiLite	Heraeus Kulzer	90 s	90 s	180 s

**Safety data sheets can be found on our website at [www.candulor.com](http://www.candulor.com).**

**PROCESSING (STEP BY STEP)****1. Preparation of tooth surface to be characterized**

Depending on the desired result to be achieved, the surface of the tooth to be characterized must be ground using a tungsten carbide bur. Attention should be paid to removing sufficient substance to allow for adequate material layers in subsequent processing steps.

**2. Conditioning of veneering surface**

At 2 bar pressure and 80–100 µm Al<sub>2</sub>O<sub>3</sub> blasting material.

**3. Cleaning**

Do **not** touch the blasted surface

- with fingers anymore
- rinse with water
- clean with steam

The pre-treated teeth are adequately wetted with monomer for 2–3 min (i.e. AESTHETIC BLUE by CANDULOR).

**4. Bond**

Stains for Resin Teeth composite is available to create the bonding layer between the resin tooth and the veneering layer to be newly characterized. The Composive is applied directly onto the surface to be characterized and distributed evenly using a brush. Adhere to layer thicknesses of min. 0.2 mm and max. 0.5 mm.

*Polymerization*

The parameters for Composive polymerization can be taken from the polymerization table.

**5. Individual characterization**

Mix the intense color stains as required and/or modify with Clear. Apply the characterizations to the polymerized Composive to suit individual requirements.

*Interim polymerization*

Polymerize finished characterized surfaces according to the polymerization table, per surface and/or quadrant.

*Attention*

Care should be taken that the shadow zones are also fully polymerized.

**6. Finisher layer**

A coating of Finisher is applied between 0.1 mm to max. 2 mm to protect the characterized surface or in case of any shape modifications.

*Polymerization*

Polymerize the modelled Finisher layer according to the parameters given in the polymerization table. Always work with an oxygen barrier layer which prevents the formation of an inhibition layer during final polymerization of the Finisher paste. The following product is suitable for this purpose: SR Gel (Ivoclar Vivadent) as well as other commercially available products. For exact application please see the instructions for use of the respective manufacturers.

**7. Finishing the polymerized surface**

Perform any shape corrections using a tungsten carbide bur. A soft rubber polisher is then suitable for pre-polishing. A further polishing step is performed using a goat hairbrush and the KMG polishing paste. Finally, carefully polish the surface with a cotton or leather buffing wheel and KGM polishing paste to achieve a high-gloss finish.

**ATTENTION**

- > In case of insufficient polish, micro-roughness may encourage deposits on the surface.
- > Improper polymerization may also lead to clinical failure.

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**CANDULOR AG**