


























## Explanation of symbols

Symbol	Symbol Titel/Description	Symbol Source <sup>[1]</sup> / Reference Number
	Manufacturer	ISO 15223 – 1 / 5.1.1
	Authorized Representative in the European Community/European Union	ISO 15223 – 1 / 5.1.2
	Date of manufacture	ISO 15223 – 1 / 5.1.3
	Use-by date	ISO 15223 – 1 / 5.1.4
	Batch code	ISO 15223 – 1 / 5.1.5
	Catalogue number	ISO 15223 – 1 / 5.1.6
	Serial number	ISO 15223 – 1 / 5.1.7
	Importeur	ISO 15223-1 / 5.1.8
	Distributor	ISO 15223 - 1 / 5.1.9
	Fragile, handle with care	ISO 15223 – 1 / 5.3.1
	Keep away from sunlight	ISO 15223 – 1 / 5.3.2
	Keep dry	ISO 15223 – 1 / 5.3.4
	Temperature limit	ISO 15223 – 1 / 5.3.7
	Do not re-use	ISO 15223 – 1 / 5.4.2



	Consult instructions for use or consult electronic instructions for use	ISO 15223 – 1 / 5.4.3
	Caution	ISO 15223 – 1 / 5.4.4
	Contains nano materials	ISO 15223 – 1 / 5.4.11
	Medical device	ISO 15223 – 1 / 5.7.7
	Unique device identifier	ISO 15223 – 1 / 5.7.10
	Health Industry Bar Code	Health Industry Business Communications Council
	CE marking is the manufacturer's declaration that the product meets the requirements of the applicable EC legislation. Where applicable: The identification number of the Notified Body shall follow this symbol.	MDR (EU) 2017/745
	To indicate that separate collection for waste electric and electronic equipment (WEEE) is required.	Directive 2012/19/EU
<b>Rx ONLY</b>	U.S. federal law restricts this device to sale by or on the order of a dental professional.	21 CFR 801.109(b) (1)
	Sterilizable in a steam sterilizer (autoclave) at the temperature specified	ISO 21531 / 20
	This way up	ISO 780 / 13
	Stacking limit by mass	ISO 780 / 15

[1]

- ISO 15223-1 Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements
- ISO 21531 Dentistry - Graphical symbols for dental instruments
- ISO 780 Packaging - Distribution packaging - Graphical symbols for handling and storage of packages